

The ibidi product family is comprised of a variety of μ-Slides and μ-Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength.

The μ-Slide I Luer I Luer is designed for cell culture under perfusion and all flow applications. Main applications are the simulation of blood vessels for arteriosclerosis research and applying defined shear stress and shear rates on cells inside the channel. The female Luers allow easy connections to tubing and pump systems. The μ-Slide I Luer comes in five versions which only differ in their channels' heights and channel volumes.

Material

ibidi μ-Slides, μ-Dishes, and μ-Plates are made of a plastic that has the highest optical quality. The polymer coverslip on the bottom exhibits extremely low birefringence and autofluorescence, similar to that of glass. Also, it is not possible to detach the bottom from the upper part. The μ-Slides, μ-Dishes, and μ-Plates are not autoclavable, since they are only temperature-stable up to 80°C/175°F. Please note that gas exchange between the medium and incubator's atmosphere occurs partially through the polymer coverslip, which should not be covered.

Optical Properties ibidi Polymer Coverslip

Refractive index n_D (589 nm)	1.52
Abbe number	56
Thickness	No. 1.5 (180 μm)
Material	polymer coverslip

Please note! The ibidi polymer coverslip is compatible with certain types of immersion oil only. A list of suitable oils can be found on page 4.

Shipping and Storage

The μ-Slides, μ-Dishes and μ-Plates are sterilized and welded in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is listed in the following table.

Conditions

Shipping conditions	Ambient
Storage conditions	RT (15-25°C)

Shelf Life of Different Surfaces

ibiTreat, Glass Bottom, ESS	36 months
Collagen, Poly-L-Lysine	18 months

Geometry of the μ-Slide I Luer

The μ-Slide I Luer provides standard slide format according to ISO 8037/1.

General Dimensions

Number of Channels	1
Channel length	50 mm
Channel width	5.0 mm
Volume per reservoir	60 μl
Growth area	2.5 cm ² per channel
Bottom matches coverslip	No. 1.5

The channel volume differs, depending on the channel height (see table below).

Product name	Channel height	Channel volume
μ-Slide I ^{0.2} Luer	200 μm	50 μl
μ-Slide I ^{0.4} Luer	400 μm	100 μl
μ-Slide I ^{0.6} Luer	600 μm	150 μl
μ-Slide I ^{0.8} Luer	800 μm	200 μl

μ-Slide Surfaces

Depending on the type of cells and the special application you are using, you will need μ-Slides with different surfaces. If you do not require any special adhesion molecules for your application, the best choice will be ibiTreat, a tissue culture treated surface.

The uncoated μ-Slide is manufactured from hydrophobic plastic. For the cultivation of most cell lines, it is indispensable to treat the uncoated μ-Slide with biopolymers, which mediate cell adhesion and growth.

The μ-Slide I Luer is also provided with a Collagen coated surface. Such an adhesion substrate has been shown to stimulate the adhesion and growth of various cell lines in μ-Slides. A high quality Collagen IV solution (Corning #356233) is used to pre-coat the slides.

Coating your μ-Slide I Luer

The uncoated μ-Slide must be coated to promote cell adhesion. If you want to establish a certain coating to match your needs, we recommend testing your coating procedure on both uncoated and ibiTreat μ-Slides, since we have observed that some biomolecules adhere differently to hydrophobic and hydrophilic plastic surfaces.

- Prepare your coating solution according to the manufacturer's specifications or reference.
- Apply the channel volume depending on the channel height (see table below) and leave at room temperature for at least 30 minutes.
- Aspirate the solution and wash with the recommended protein dilution buffer.
- Optionally let dry at room temperature. Attention, some coating proteins might degenerate when drying!

Detailed information about coatings is provided in [Application Note 08 Cell culture coating](#).

Product Name	Channel Volume	Coating Area
μ-Slide I ^{0.2} Luer	50 μl	5.2 cm ²
μ-Slide I ^{0.4} Luer	100 μl	5.4 cm ²
μ-Slide I ^{0.6} Luer	150 μl	5.6 cm ²
μ-Slide I ^{0.8} Luer	200 μl	5.8 cm ²


Tip:

You can add the buffer into one channel end and simultaneously aspirate it on the other side.


Filling and Handling of Channel Slides

In order to avoid air bubbles inside the channels please follow the recommendations below.

When filling the channels put the pipet tip directly to the channel's inlet. Apply the volume with a constant and swift flow.




In special cases, e.g. when the channel surface is hydrophobic or when filling small channels, it might be necessary to fill the channel with a syringe. Use a low volume syringe with 1 or 2.5 ml!



Important!

When seeding cells, fill only the correct channel volume into the channel. Avoid surplus cell suspension in the reservoirs!


After cell attachment fill 60 μl in each well, for a better medium supply to the cells.



Medium Exchange

The following medium exchange protocol is important for cell culture medium exchange, staining, washing and coating procedures.


Empty the reservoirs completely without emptying the channel. Inject the new solution from one side and remove the old solution from the other side. Make sure the old solution is completely replaced. For a 99% exchange add about three times the channel volume from one side.



Important!

Take care, that the channel never falls dry during the exchange process. This helps you avoiding air bubbles!

When aspirating the liquid put the pipet tip away from the channel's inlet! This prevents you from evacuating the whole channel.



Cell Culture under Static Conditions

For many static applications with microscopic imaging, like transfection, immunofluorescence staining or cell morphology the μ-Slide I Luer is an optimal solution.

Important!

The μ-Slide I^{0.2} Luer is not recommended for use in static cell culture!
For longer cultivation, a gentle flow is necessary. This can be achieved by a perfusion system or an incubator-compatible cell culture rocker.

- Trypsinize and count cells as usual. The cell density after seeding strongly depends on the channel's height. We recommend the following cell concentrations and volumes:

Product name	Volume	Cell concentration
μ-Slide I ^{0.2} Luer	50 μl	6–14 x 10 ⁵ cells/ml
μ-Slide I ^{0.4} Luer	100 μl	3–7 x 10 ⁵ cells/ml
μ-Slide I ^{0.6} Luer	150 μl	2–4.5 x 10 ⁵ cells/ml
μ-Slide I ^{0.8} Luer	200 μl	1.5–3.5 x 10 ⁵ cells/ml

- Apply the volume directly into the channel. The recommended cell concentration should result in a 50 % optical confluence layer after 24 hours.
- Cover reservoirs with the supplied caps. Incubate at 37°C and 5% CO₂ as usual.
- After cell attachment fill each reservoir with 60 μl medium.

Depending on the cells we recommend exchanging the medium every day in static culture: Aspirate both reservoirs (not the channel). Flush fresh medium inside the channel by filling one reservoir with 120 μl medium and removing the content of the reservoir from the other well, ensuring the channel is never dry. Leave both reservoirs filled with approx. 60 μl each.

Tip:

The day before seeding the cells we recommend placing the cell medium, the μ-Slide, and the tubing into the incubator for equilibration. This will prevent the liquid inside the channel from emerging air bubbles over the incubation time.
Quick dispensing of cell suspension helps to avoid trapped air bubbles and leads to maximal homogeneity of cell distribution.

Cell Culture under Flow Conditions

Due to the Luer adapters, μ-Slide I Luer is suitable to any fluidic setup for cell cultivation under flow conditions. Cells are seeded into the channel and the flow is applied after cell attachment.

- Trypsinize and count cells as usual. The cell density after seeding strongly depends on the channel's height. We recommend the following cell concentrations and volumes:

Product name	Volume	Cell concentration
μ-Slide I ^{0.2} Luer	50 μl	2.5–5 × 10 ⁶ cells/ml
μ-Slide I ^{0.4} Luer	100 μl	1.2–2.5 × 10 ⁶ cells/ml
μ-Slide I ^{0.6} Luer	150 μl	0.8–1.6 × 10 ⁶ cells/ml
μ-Slide I ^{0.8} Luer	200 μl	0.6–1.2 × 10 ⁶ cells/ml

- Apply the volume directly into the channel. The recommended cell concentration should result in a 100 % optical confluence layer after some hours.
- Cover reservoirs with the supplied caps. Incubate at 37°C and 5 % CO₂ as usual.
- After cell attachment fill each reservoir with 60 μl medium.
- The μ-Slide is now ready for applying flow conditions on the adherent cells. Don't trap air bubbles when plugging in the connecting tubes.

Tip:

The day before seeding the cells we recommend placing the cell medium, the μ-Slide, and the tubing into the incubator for equilibration. This will prevent the liquid inside the channel from emerging air bubbles over the incubation time.

Quick dispensing of cell suspension helps to avoid trapped air bubbles and leads to maximal homogeneity of cell distribution.

For long term analysis of cells under flow conditions we recommend using μ-Slides with ibiTreat surface.


[Application Note 13 "HUVECs under perfusion"](#) describes a detailed protocol of a long term experiment with HUVECs and the ibidi Pump System.

Detailed information about flow rates, shear stress, and shear rates is provided in [Application Note 11 "Shear stress and shear rates"](#) on www.ibidi.com.

For connecting several μ-Slides I Luer with each other in a serial way, please refer to our [Application Note 25 "Serial](#)

Connection of Flow Chamber".

Suitable Tube Adapter Sets are also available (see page 5). They consist of a tubing (20 cm) with inner diameter of 1.6 mm and adapters for the connection between the ibidi μ-Slide (female Luer) and the tubing of the pump in use.



Please contact us for recommended perfusion setups. ibidi provides a variety of channel slides and pump systems.

Preparation for Cell Microscopy

To analyze your cells, no special preparations are necessary. Cells can be observed live, or fixed directly in the μ-Slide on an inverted microscope. You can use any fixative of your choice. The μ-Slide material is compatible with a variety of chemicals, e.g., acetone or methanol. Further specifications can be found at www.ibidi.com. Due to the thin bottom of only 180 μm, high resolution microscopy is possible.


Immersion Oil

When using oil immersion objectives, use only the immersion oils specified in the table. The use of a non-recommended oil could lead to the damage of the plastic material and the objective.


Company	Product	Ordering Number
Zeiss	Immersol 518 F	(Zeiss) 444960
Zeiss	Immersol W 2010	(Zeiss) 444969
Leica	Immersion liquid	(Leica) 11513859

Ordering Information


The μ-Slide I Luer is available in different channel heights and surfaces. See table below for choosing your μ-Slide I Luer.

μ-Slide I^{0.2} Luer



Cat. No.	Description
80166	μ-Slide I ^{0.2} Luer ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
81162	μ-Slide I ^{0.2} Luer Collagen IV: #1.5 polymer coverslip, sterilized
81161	μ-Slide I ^{0.2} Luer Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized

μ-Slide I^{0.4} Luer



Cat. No.	Description
80176	μ-Slide I ^{0.4} Luer ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
81172	μ-Slide I ^{0.4} Luer Collagen IV: #1.5 polymer coverslip, sterilized
81171	μ-Slide I ^{0.4} Luer Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized

μ-Slide I^{0.6} Luer


Cat. No.	Description
80186	μ-Slide I ^{0.6} Luer ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
81182	μ-Slide I ^{0.6} Luer Collagen IV: #1.5 polymer coverslip, sterilized
81181	μ-Slide I ^{0.6} Luer Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized

μ-Slide I^{0.8} Luer


Cat. No.	Description
80196	μ-Slide I ^{0.8} Luer ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
81192	μ-Slide I ^{0.8} Luer Collagen IV: #1.5 polymer coverslip, sterilized
81191	μ-Slide I ^{0.8} Luer Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized

Tube Adapter Set


Cat. No.	Description
10831	Tube Adapter Set: sterilized

Selected References

- M. Anton, A. Wolf, O. Mykhaylyk, C. Koch, B. Gansbacher, and C. Plank. Optimizing Adenoviral transduction of endothelial cells under flow conditions. *Pharmaceutical Research*, 2012. doi: 10.1007/s11095-011-0631-2.
- Pablo Bartolucci, Vicky Chaar, Julien Picot, Dora Bachir, Anoosha Habibi, Christine Fauroux, Frederic Galacteros, Yves Colin, Caroline Le Van Kim, and Wassim El Nemer. Decreased sickle red blood cell adhesion to laminin by hydroxyurea is associated with inhibition of lu/bcam protein phosphorylation. *Blood*, 2010. doi: 10.1182/blood-2009-12-257444.
- C. Boehlke, F. Kotsis, V. Patel, S. Braeg, H. Voelker, S. Bredt, T. Beyer, H. Janusch, C. Hamann, M. Gödel, K. Müller, M. Herbst, M. Hornung, M. Doerken, M. Köttgen, R. Nitschke, P. Igarashi, G. Walz, and E. W. Kuehn. Primary Cilia Regulate mTORC1 Activity and Cell Size Through Lkb1. *Nature Cell Biology*, 2010. doi: 10.1038/ncb2117.
- A. Sabine, Y. Agalarov, H. Maby-El Hajjami, M. Jaquet, R. Hägerling, C. Pollmann, D. Bebbber, A. Pfenniger, N. Miura, O. Dormond, J.-M. Calmes, Ralf H. Adams, T. Mäkinen, F. Kiefer, B. Kwak, and T. Petrova. Mechanotransduction, PROX1, and FOXC2 Cooperate to Control Connexin37 and Calcineurin during Lymphatic-Valve Formation. *Developmental Cell*, 2012.
- T. G. Walsh, R. P. Murphy, P. Fitzpatrick, K. D. Rochfort, A. F. Guinan, A. Murphy, and P. M. Cummins. Stabilization of Brain Microvascular Endothelial Barrier Function by Shear Stress Involves VE cadherin Signaling Leading to Modulation of pTyr Occludin Levels. *Journal of Cellular Physiology*, 2011. doi: 10.1002/jcp.22655.

For research use only!

Further technical specifications can be found at www.ibidi.com. For questions and suggestions please contact us by e-mail info@ibidi.de or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany.

© ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.