

The ibidi product family is comprised of a variety of μ-Slides and μ-Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength.

The μ-Slide y-shaped can be easily connected to a pump via its Luer connectors. It enables you to grow cells under flow conditions with a bifurcation of 30 and 45 degree depending on the flow direction you choose. It is meant as a simulation system for blood vessels where the reaction of cells to a stimulus of your choice can be observed in real-time.

## Material

ibidi μ-Slides, μ-Dishes, and μ-Plates are made of a plastic that has the highest optical quality. The polymer coverslip on the bottom exhibits extremely low birefringence and autofluorescence, similar to that of glass. Also, it is not possible to detach the bottom from the upper part. The μ-Slides, μ-Dishes, and μ-Plates are not autoclavable, since they are only temperature-stable up to 80°C/175°F. Please note that gas exchange between the medium and incubator's atmosphere occurs partially through the polymer coverslip, which should not be covered.

### Optical Properties ibidi Polymer Coverslip

Refractive index $n_D$ (589 nm)	1.52
Abbe number	56
Thickness	No. 1.5 (180 μm)
Material	polymer coverslip

**Please note! The ibidi polymer coverslip is compatible with certain types of immersion oil only. A list of suitable oils can be found on page 3.**

## Shipping and Storage

The μ-Slides, μ-Dishes and μ-Plates are sterilized and welded in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is listed in the following table.

### Conditions

Shipping conditions	Ambient
Storage conditions	RT (15-25°C)

### Shelf Life of Different Surfaces

ibiTreat, Glass Bottom, ESS	36 months
Collagen, Poly-L-Lysine	18 months

## Geometry of the μ-Slide y-shaped

The μ-Slide y-shaped provides a standard slide format according to ISO 8037/1.

### Dimensions

Channel volume	110 μl
Channel height	0.4 mm
Branching angles	30° and 45°
Adapters	female Luer
Volume per reservoir	60 μl
Growth area	2.8 cm <sup>2</sup>
Coating area using 110 μl	5.6 cm <sup>2</sup>
Bottom matches coverslip	No. 1.5

## μ-Slide Surfaces

Depending on the type of cells and the special application you are using, you will need μ-Slides with different surfaces. If you do not require any special adhesion molecules for your application, the best choice will be ibiTreat, a tissue culture treated surface.

The uncoated μ-Slide is manufactured from hydrophobic plastic. For the cultivation of most cell lines, it is indispensable to treat the uncoated μ-Slide with biopolymers, which mediate cell adhesion and growth.

## Coating your μ-Slide y-shaped

The uncoated μ-Slide must be coated to promote cell adhesion. If you want to establish a certain coating to match your needs, we recommend testing your coating procedure on both uncoated and ibiTreat μ-Slides, since we have observed that some biomolecules adhere differently to hydrophobic and hydrophilic plastic surfaces.

- Prepare your coating solution according to the manufacturer's specifications or reference.

- Apply 110 μl and leave at room temperature for at least 30 minutes.
- Aspirate the solution and wash with the recommended protein dilution buffer.
- Optionally let dry at room temperature. Attention, some coating proteins might degenerate when drying!

Detailed information about coatings is provided in [Application Note 08 Cell culture coating](#).

**Tip:**

Apply the coating solution with a 1000 μl pipette. Adding more liquid than needed (e.g. 200 μl) might be more convenient. Aspirate the surplus 90 μl from the Luer reservoirs after filling.  
For washing you can add the buffer into one channel end and simultaneously aspirate it on the other side.

**Tip:**

The day before seeding the cells we recommend to place the cell medium and the μ-Slide into the incubator for equilibration. This will prevent the liquid inside the channel from emerging air bubbles over the incubation time.

Quick dispensing of cell suspension helps to avoid trapped air bubbles and leads to maximal homogeneity of cell distribution.

For long term analysis of cells under flow conditions we recommend to use μ-Slides with the ibiTreat surface.

**Seeding Cells**

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, application of a  $3-7 \times 10^5$  cells/ml suspension should result in a confluent layer within 2-3 days.
- Apply 200 μl cell suspension into the channel of the μ-Slide using a 1000 μl pipette.
- Remove the cell suspension from the reservoirs.
- Cover reservoirs with the supplied caps. Incubate at 37°C and 5% CO<sub>2</sub> as usual.
- After cell attachment fill each reservoir with 60 μl medium for longer cultivation.
- The μ-Slide is now ready for applying flow conditions on the adherent cells. Don't trap air bubbles when plugging in the connecting tubes.

Depending on the cells we recommend to exchange the medium every day in static culture: Aspirate both reservoirs. Flush fresh medium inside the channel by slowly filling one reservoir with 400 μl medium and removing the content of the reservoir from the other well, ensuring the channel is never dry. Leave both reservoirs filled with approx. 60 μl each.


**Preparation for Cell Microscopy**

To analyze your cells, no special preparations are necessary. Cells can be observed live, or fixed directly in the μ-Slide on an inverted microscope. You can use any fixative of your choice. The μ-Slide material is compatible with a variety of chemicals, e.g., acetone or methanol. Further specifications can be found at [www.ibidi.com](http://www.ibidi.com). Due to the thin bottom of only 180 μm, high resolution microscopy is possible.

**Flow Application**

Detailed information about flow rates, shear stress, and shear rates is provided in [Application Note 11 "Shear stress and shear rates"](#) and [Application Note 18 "Shear Stress and Shear Rates in μ-Slide y-shaped"](#).

Suitable Tube Adapter Sets are also available (see page 4). They consist of a tubing (20 cm) with inner diameter of 1.6 mm and adapters for the connection between the ibidi μ-Slide (female Luer) and the tubing of the pump in use.



Please contact us for recommended perfusion setups. ibidi provides a variety of channel slides and pump systems.


**Immersion Oil**

When using oil immersion objectives, use only the immersion oils specified in the table. The use of a non-recommended oil could lead to the damage of the plastic material and the objective.

Company	Product	Ordering Number
Zeiss	Immersol 518 F	(Zeiss) 444960
Zeiss	Immersol W 2010	(Zeiss) 444969
Leica	Immersion liquid	(Leica) 11513859


**Ordering Information**

The μ-Slide y-shaped is available with different surfaces. See table below for choosing your μ-Slide y-shaped.



Cat. No.	Description
80126	μ-Slide y-shaped ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
80121	μ-Slide y-shaped Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized

**Tube Adapter Set**



Cat. No.	Description
10831	Tube Adapter Set: sterilized

**Selected References**

- N. Cockcroft, O. Oke, F. Cunningham, E. Bishop, I. M. Fearon, R. Zantl, and M. D. Gaça. An In Vitro Perfusion System to Examine the Responses of Endothelial Cells to Simulated Flow and Inflammatory Stimulation. *ATLA*, 2009.
- J. Samarin, I. Cicha, and M. Goppelt-Struebe. Cell type-specific regulation of CCN2 protein expression by PI3K-AKT-FoxO signaling. *Journal of Cell Communication and Signaling*, 2009. doi: 10.1007/s12079-009-0055-5.
- M. J. Seidler, S. Salvenmoser, and F.-M. C. Muller. *Aspergillus fumigatus* Forms Biofilms with Reduced Antifungal Drug Susceptibility on Bronchial Epithelia Cells. *Antimicrob. Agents Chemother.*, 2008. doi: 10.1128/aac.00234-08.
- K. Urschel, C. D. Garlichs, W. G. Daniel, and I. Cicha. VEGFR2 signalling contributes to increased endothelial susceptibility to TNF-α under chronic non-uniform shear stress. *Atherosclerosis*, 2011.

**For research use only!**

Further technical specifications can be found at [www.ibidi.com](http://www.ibidi.com). For questions and suggestions please contact us by e-mail [info@ibidi.de](mailto:info@ibidi.de) or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany.

© ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.