

# SCANSPEED 406

## TECHNICAL DATA



### HIGHLIGHTS

- Compact & robust small centrifuge, which is ideal for routine low-speed centrifugation.
- Commonly used in the fractionation of blood samples & precipitation of biological samples.
- Automatic lid release and vivid end alarm for easy retrieval of samples.
- Supplied with an angle rotor including 6 aluminium sleeves for 15 mL tubes.

For more information please visit  
[www.labogene.com](http://www.labogene.com)

### Details

Product family	ScanSpeed
Model:	406
Part number	7607514101
Product type	Clinical low-speed centrifuge
Compliance with	EN 61010-1, EN 61010-2-020, EN 61000-6-1, EN 11201
e-catalogue classification (UNSPSC code)	40161701 Centrifuges

### Technical specifications

### Unit

Typical application		Routine separation of blood, urine or cells.
Max. RPM/RCF	RPM/x g	4.000/2.075
Max capacity		6 x 15 mL
Time control		Time < 100 min or continuous
RPM/RCF conversion		Yes
Noise level	dB(A)	≤52
ACC/DEC	sec	≤20 / ≤20
Program memory		10
Parameters on display window		RPM(RCF) & min:sec
Imbalance cutout		Yes
Safety lid lock		Yes
Lid drop protection		Yes
Automatic lid release at completion		Yes
Power supply	V/Hz	220/50-60
Power requirements	VA	140
Dimensions (WxDxH)	mm	296x412x206
Weight without rotor	kg	17,5

**VITARIS**  
 A Difference in Life Science

**VITARIS AG** Tel.: 041 769 00 00  
 Blegistrasse 9 Fax: 041 769 00 01  
 CH-6340 Baar E-mail: [info@vitaris.com](mailto:info@vitaris.com)

[www.vitaris.com](http://www.vitaris.com)

**LABOGENE**  
 Scandinavian by Design