

Is your safety cabinet still in shape?

Don't leave anything to chance.

With our FITCheck you have safety in black and white.

Regular inspection of safety cabinets is mandatory for compliance with legal regulations. A valid air flow and intact filters are essential to ensure both personal protection and product integrity, and to detect malfunctions at an early stage.

In addition, regular tests can reduce energy costs and increase comfort.



995.-

excl. quantity rate

Scope of services of the FITCheck:

- In-depth inspection by certified service technicians
- Functional control and adjustment of all safety, control and operating elements as well as mechanical and electrical components
- Air velocity measurement of laminar airflow, downflow and inflow
- Leak testing of main and exhaust air HEPA filters according to EN 12469
- ISO 5 cleanliness class test
- Maintenance log with printout of all measured values – clear and documented.

Your benefits:

- **Regulatory Compliance**
- **Protection against unforeseeable failures**
- **Energy efficiency and cost reductions**
- **Preservation of the protective effect**
- **Enhanced user experience**
- **If a new safety cabinet is purchased as a result of the FITCheck, the check price will be offset against the purchase price**

FITCheck offer: CHF 995.- per device

flat rate incl. VAT

Volume discounts for the same location:

- 2 devices: 5% discount
- 3 devices: 10% discount
- 4 devices: 15% discount
- from 5 devices: 20% discount

Call and make an appointment: 041 769 00 75