

Cool Deal: PHCBI MPR-715PF-PE

phcbi

Only a few units left!

The MPR-715F-PE combines high-performance refrigeration with a -30°C freezer compartment in a single unit.

With 415 liters of refrigerated capacity and 176 liters of freezer capacity, the MPR-715F-PE offers flexible, reliable, and convenient sample storage in a compact footprint.



CHF 6,110
instead of CHF 10,590
while supplies last

- **Optimal Temperature Distribution**
The air circulation generated by a fan in the refrigerator compartment ensures excellent uniformity and rapid temperature recovery after routine door openings. This helps minimize fluctuations and protects temperature-sensitive materials.
- **Precise Control**
A microprocessor and precise sensors maintain accurate temperatures in both the refrigerator and freezer compartments. A comprehensive alarm system alerts users to temperature deviations, door openings, and power outages.
- **User-Friendliness Meets Performance**
The four-door design reduces air loss when the door is opened and makes efficient use of the space in front of the unit. LED lighting, adjustable shelves, and key locks support practical daily use in demanding laboratory environments.

Take advantage now: info@vitaris.com